# Female Ultrasound Training Pelvis

#### Model 404A



# **TEACH • TRAIN • PRACTICE**

The gynecological ultrasound exam is an important diagnostic tool. Students gain competency with hands-on practice. Unfortunately, access to patients can be limited. The Model 404A Female Ultrasound Training Pelvis facilitates teaching and demonstration of gynecological ultrasound techniques in a non-stressful situation.

Our Model 404A consists of an external female pelvic model containing a uterus, fallopian tubes, ovaries, bladder and rectal landmarks useful for scanning. The phantom has both abdominal and vaginal scanning access allowing a variety of transducer orientations.

The Female Ultrasound Training Pelvis creates a relaxed learning environment for teaching and developing ultrasound examination skills and techniques as well as demonstrating 3D ultrasound capabilities.

#### Features

- External anatomy
- Abdominal and vaginal scanning access
- Uterus with endometrium and myometrium
- Ovary with follicles
- Full bladder
- · Rectal landmarks
- Realistic Scanning



www.3000buy.com 深圳市净康科技有限公司

### FEMALE ULTRASOUND TRAINING PELVIS

## Model 404A



**Transabdominal of uterus** 

## **SPECIFICATIONS**



Trans-vaginal of uterus, bladder and rectum



**Ovary with follicles** 

| DIMENSIONS         | 15 cm x 11 cm x 11 cm                                       |
|--------------------|---|
|                    | (5.9" x 4.3" x 4.3")  |
| PHANTOM WEIGHT     | 2 kg (4.4 lbs.)   |
| MATERIALS          | Z-Skin <sup>™</sup> elastomer and Zerdine®                  |
| INTERNAL STRUCTURE | Vaginal Canal   |
|                    | Cervix  |
|                    | Uterus with Endometrium, Myome-<br>trium and Uterine Cavity |
|                    | Fallopian Tubes   |
|                    | Homogeneous Left Ovary                                      |
|                    | Right Ovary containing follicles                            |
|                    | Iliac Vessels   |

#### **MODEL 404A INCLUDES**

| QTY   | COMPONENT DESCRIPTION             |
|---|-----------------------------------|
| 1   | Female Ultrasound Training Pelvis |
| 1   | User Guide                        |
| -   | 12-Month Warranty*                |
| *Please note: once any device has been inserted into the phantom (such as |                                   |

biopsy needles, localization wires, etc.) the 12 month warranty will not cover claims related to material desiccation or needle tracking.

